

Product Change Notification

Subject: Circuit change for Dino-Lite® Premier AM4113 and AD4113 series

Products Affected:

All Dino-Lite Premier AM4113 and AD4113 series

Description of Change:

The circuit, including sensor, will be changed to improve reliability and image quality of these series. There will be no difference on the product name, physical appearance, accessories, packing, and price. For identification purposes, the product name on the label will be noted with “(R4)” which means Release 2014, e.g. AD4113 (R4). The product name on other information sources, such as the website, will not change.

For SDK user’s attention, the AM4113 (R4) and AD4113 (R4) add MJPG mode. Although the 4113 and 4113 (R4) series models share the same SDK commands, we advise paying special attention to the color and exposure conditions, as they may be affected.

Effective Date

The AM4113 (R4) and AD4113 (R4) series will be available starting from June 2014.



For any questions regarding this notification, please contact your AnMo contact.